


Technical Specifications

| Commercial name | MICRO F | | POD F | |
|-----------------------------------|---|----------------|--|----------------|
| Material | 25% hydrophilic acrylic | | 26% hydrophilic acrylic | |
| Overall diameter | 10.75 mm | | 11.40 mm | |
| Optic diameter | 6.15 mm | | 6.00 mm | |
| Optic | Biconvex aspheric (-0.11 μ SA) trifocal diffractive FINEVISION | | | |
| Filtration | UV & blue light | | | |
| Refractive index | 1.46 | | | |
| Abbe number | 58 | | | |
| Angulation | 5° | | | |
| Additional power | + 1.75D for intermediate vision and + 3.50D for near vision | | | |
| Injection system | Medicel Viscoject Bio 1.8/2.2 up to 24.5D Medicel Accuject 1.8/2.0 up to 24.5D Medicel Accuject 2.1/2.2 up to 35D | | Medicel Accuject 2.0 up to 24.5D Medicel Accuject 2.1/2.2 up to 35D | |
| Incision size | ≥ 1.8 mm | | ≥ 2.0 mm | |
| Spherical power | 10D to 35D (0.5D steps) | | 6D to 35D (0.5D steps) | |
| Square edge | 360° | | | |
| Nominal manufacturer A constant | 118.80 | | 118.95 | |
| Suggested A constant ¹ | Interferometry | Ultrasound | Interferometry | Ultrasound |
| Hoffer Q: pACD | 5.35 | 5.26 | 5.59 | 5.35 |
| Holladay 1: Sf | 1.60 | 1.48 | 1.83 | 1.57 |
| Barrett: LF | 1.78 | - | 1.86 | - |
| SRK/T: A | 118.80 | 118.59 | 118.95 | 118.73 |
| Haigis ² : a0; a1; a2 | 1.36; 0.4; 0.1 | 1.13; 0.4; 0.1 | 1.36; 0.4; 0.1 | 1.13; 0.4; 0.1 |

¹ Estimates only; surgeons are recommended to use their own values based upon their personal experience. Refer to our website for updates.

² Not optimized.

Product Information

| | MICRO F | POD F |
|--------------------------------|--|--|
| Manufacturer | PhysIOL s.a. - Liège Science Park Allée des Noisetiers 4 B-4031 Belgium +32 4 361 05 49 physiol@bvimedical.com | PhysIOL s.a. - Liège Science Park Allée des Noisetiers 4 B-4031 Belgium +32 4 361 05 49 physiol@bvimedical.com |
| Certificate information | CE: Certificate N° CE658516 ISO 13485:2016: Certificate n° MD658518 MDSAP: Certificate N° MDSAP 691544 ISO 9001:2015: Certificate N° FM 658519 | CE: Certificate N° CE658516 ISO 13485:2016: Certificate n° MD658518 MDSAP: Certificate N° MDSAP 691544 ISO 9001:2015: Certificate N° FM 658519 |
| Shelf life | Five (5) years from manufacturing date | Five (5) years from manufacturing date |
| Intended Use | Intended use (for all IOLs): The posterior chamber intraocular lens which is intended to be placed into the capsular bag for the replacement of the human lens to achieve the visual correction of aphakia in adult patients in whom the cataractous lens has been removed by extracapsular cataract extraction. | Intended use (for all IOLs): The posterior chamber intraocular lens which is intended to be placed into the capsular bag for the replacement of the human lens to achieve the visual correction of aphakia in adult patients in whom the cataractous lens has been removed by extracapsular cataract extraction. |
| Indication for use | The lens should be used as intended in patients surgically treated for cataract, with possibly associated presbyopia, who desire improved uncorrected far vision, useful near and intermediate visual functions and reduced spectacle dependence. | The lens should be used as intended in patients surgically treated for cataract, with possibly associated presbyopia, who desire improved uncorrected far vision, useful near and intermediate visual functions and reduced spectacle dependence. |
| Product Composition | No products of animal or human origin are present in the implant. The implant is made of the HELIO25 material, composed of an acrylate copolymer Hydroxyethyl methacrylate (HEMA) and Ethoxyethyl methacrylate (EOEMA), including a UV and blue light filter | No products of animal or human origin are present in the implant. The implant is made of the HELIOFLEX material, composed of an acrylate copolymer Hydroxyethyl methacrylate (HEMA) and Methyl methacrylate (MMA), including a UV and blue light filter |
| For sterile product | All IOLs from PhysIOL are steam sterilized | All IOLs from PhysIOL are steam sterilized |
| Packaging Material | Holder (Polypropylene) Container (Polypropylene) Storage liquid (0.9% NaCl solution) Aluminium lid (Aluminium Gold) Container label (paper) Blister PP (Polypropylene) Tyvek lid | Holder (Polypropylene) Container (Polypropylene) Storage liquid (0.9% NaCl solution) Aluminium lid (Aluminium Gold) Container label (paper) Blister PP (Polypropylene) Tyvek lid |
| Product Class | MDD Class IIb Sterile, According to European Medical Device Directive 93/42/EEC Not available in the United States | MDD Class IIb Sterile, According to European Medical Device Directive 93/42/EEC Not available in the United States |

Injection Guidelines

The Medical Viscoject Bio 1.8/2.2 and Accuject 1.8/2.0/2.1/2.2 injection system are recommended for implanting the FINEVISION lenses. This fully single-use system represents total reliability for safe and effective lens injections.


Its compact design with integrated cartridge enables a simple, predictable loading and positioning of the lens.

FINEVISION MICRO F

Medical Viscoject Bio 1.8/2.2 up to 24.5D

Medical Accuject 1.8/2.0 up to 24.5D and 2.1/2.2 up to 35D

Guidelines with Viscoject:




1. Apply ophthalmic viscoelastic device (OVD) into the tip and the loading chamber of the injector cartridge.
2. Remove the lens from the lens holder. Position the lens into the cartridge in such way that the two haptics with the notches are pointing at 1 and 7 o'clock.
3. Exert slight pressure onto the lens optic and make sure that all haptics are inside before further closing the cartridge.
4. Close the cartridge and check the position of the lens. Once the "click-lock" mechanism engages, the lens is securely loaded.
5. Engage the cartridge in the injector.
6. Press the injector plunger forward and push the lens into the conical tip of the cartridge. Pull the plunger back a few millimeters and then inject the lens in one continuous motion. For gentle implantation, it is not necessary to fully push the plunger to the bottom of the cartridge.

FINEVISION POD F

Accuject 2.0 for lens diopters up to 24.5D

Accuject 2.1/2.2 for lens diopters up to 35D

Guidelines with Accuject:



1. Apply ophthalmic viscoelastic device (OVD) into the tip and the loading chamber of the injector cartridge.
2. Remove the lens from the lens holder. Position the lens into the cartridge in such way that the two haptics with the notches are pointing at 1 and 7 o'clock.
3. Exert slight pressure onto the lens optic and make sure that all haptics are inside before further closing the cartridge. Close the cartridge and check the position of the lens.
4. Once the "click-lock" mechanism engages, the lens is securely loaded and ready for injection.
5. Press the injector plunger forward and push the lens into the conical tip of the cartridge.
6. Pull the plunger back a few millimeters and then inject the lens in one continuous motion. For gentle implantation, it is not necessary to fully push the plunger to the bottom of the cartridge.